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SV Series Vibro Viscometer
Water Jacket Assembly

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© 2004 A&D Company Ltd. All rights reserved.

No part of this publication may be reproduced, transmitted, transcribed, or translated into any language in any form by any means without the written permission of A&D Company Ltd.

The contents of this manual and the specifications of the instrument covered by this manual are subject to change for improvement without notice.

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ENGLISH

1. Introduction

The AX-SV-37, a water jacket assembly, is used with the SV Series Vibro Viscometer, to maintain the sample temperature constant or to measure the viscosity when the sample temperature is changed. The operating temperature range of the AX-SV-37 is 0°C to 100°C

The AX-SV-37 consists of the following:

Water jacket	1 pc	(Main body: PC (Polycarbonate), Packing: Silicone rubber, Washer: Nylon)
Small sample cup (10 mL)	4 pcs	(PC (Polycarbonate), Operating temperature: 120°C max.)
Small sample cup cover	4 pcs	(PC (Polycarbonate), Operating temperature: 120°C max.)
Screw	1 pc	(Screw: POM (Polyacetal), Washer: Nylon)

Note

- As a heating medium, use water, isopropyl alcohol (IPA), mixture of water and IPA, or silicone oil, which does not erode the materials described above. Using the heating medium other than these may damage the water jacket.
- When using the water jacket, make sure that no inner pressure is exerted in the water jacket due to the kinked or blocked tubes, as that could exert the pressure in the water jacket, causing it to break.
- Set the circulation flow rate to 5 L/min or less. Otherwise, the water jacket assembly may be damaged.
- To control the temperature, a commercially available constant temperature bath is required separately.


2. Installation

As a sample container, the small sample cup provided or the glass sample cup (AX-SV-35) sold separately as an option, is used.

To circulate the heating medium into the water jacket, a commercially available constant temperature bath is required.

1. Using silicone tubes with an inside diameter of 8 mm, make a connection between the "OUT" connector of the bath and the lower nozzle of the water jacket indicated as **1** in the illustration, and then, between the "IN" connector and the upper nozzle indicated as **2**.

Note: It is recommended that, for safety, a commercially available tubing clamp (clamping size: 11 to 20 mm) be used to fasten the tube securely to the nozzle.

2. Pour the heating medium into the water jacket indicated as **3**. An appropriate sample amount is approximately 11 mL.

The heating medium conducts the heat of the water jacket to the sample cup.


Note: When the levels of the sample surface and the heating medium surface are the same and hard to distinguish the positioning of the sensor plates, change the amount of the heating medium.

3. Pour the sample into the sample cup. When the small sample cup is used, pour the sample up to the 10 mL line.

4. Insert the sample cup into the water jacket indicated as **3**.


When the specific gravity of the sample liquid is small and the small sample cup floats, secure lightly the front side of the cup to the water jacket using the screw provided.

Use the small sample cup cover for volatile samples.


If the handles of the small sample cup interfere with the measurement, they can be removed.

Apply force in the up and down directions slowly as shown in the illustration, to break the handle off.


3. How to Use

- Attach the water jacket on the table along the guides. Make sure that the water jacket is slid furthest along the guides until stopped.


- When the water jacket is used, be sure to raise or remove the protector.

- Use much care so that the sensor plates and the temperature sensor do not come into contact with the inner wall of the sample. If they are in contact, a measurement error may occur.

- Be sure to adjust the sample surface to the center of the narrow part of the sensor plates.

4. Measuring the Absolute Value of Viscosity Using the Small Sample Cup

The viscometer has been calibrated using the accessory sample cup (35 mL) when shipped. The distance between the inner wall and the sensor plates when the accessory sample cup is used, differs from the distance when the small sample cup (10 mL) is used. This causes a difference in the sensor plate's detection capability, thus causing a difference in the viscosity measured.

Therefore, to measure the absolute value of viscosity using the small sample cup (10 mL), it is recommended that calibration be performed using a fluid with a known viscosity value which is close to the sample viscosity. (Refer to "6. VISCOSITY CALIBRATION" of the viscometer instruction manual.)

5. Maintenance


Clean the small sample cup as necessary. If cleaning is not sufficient, a measurement error may occur due to contamination.

6. Specifications

Operating temperature:	0°C to 100°C
Circulation nozzle outside diameter:	10.5 mm
Recommended hose:	Silicone tube with an inside diameter of 8 mm

- When an stirrer is attached at the bottom of the water jacket, the sample viscosity can be measured while the sample is being stirred. The maximum viscosity value which can be measured is 1000 mPa·s.

Stirrer: VARIOMAG MICRO manufactured by H+P Labortechnik AG

- Note:
- Use a rotator with a size of 6 mm (length) x 4 mm (diameter).
 - Attach the stirrer to the bottom of the water jacket, using an M4 screw from below.
 - When the rotator is used, make sure that the rotator will not touch the temperature sensor or the sensor plates. If it does, increase the sample amount so that the rotator will be farther away from the sensor plates.
 - Use the stirrer with a revolution speed which will not cause ripples in the sample surface.

THIS PAGE INTENTIONALLY LEFT BLANK.

A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013, JAPAN
Telephone: [81] (3) 5391-6132 Fax: [81] (3) 5391-1566

A&D ENGINEERING, INC.

47747 Warm Springs Blvd, Fremont, California 94539, U.S.A.
Tel: [1] (800) 726-3364 Weighing Support:[1] (888) 726-5931 Inspection Support:[1] (855) 332-8815

A&D INSTRUMENTS LIMITED

Unit 24/26 Blacklands Way, Abingdon Business Park, Abingdon, Oxfordshire OX14 1DY United Kingdom
Telephone: [44] (1235) 550420 Fax: [44] (1235) 550485

A&D AUSTRALASIA PTY LTD

32 Dew Street, Thebarton, South Australia 5031, AUSTRALIA
Telephone: [61] (8) 8301-8100 Fax: [61] (8) 8352-7409

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A&D Instruments India Private Limited

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A&D KOREA Limited