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
Instruction Manual for AD-1675A Tabletop Breeze Break (S)

1. Features


- The tabletop breeze break protects the balance from breezes due to air conditioning and body motion, thereby reducing balance weighing error.
- The transparent panel assembly consists of antistatic plastic material, which protects the balance from static electricity.
- The hole on the back panel can be used to pass the RS-232C and AC adapter cables through.
- Sliding panels on top, right and left side enable samples to be inserted or removed from each side.

2. Specifications

- Weight : Approximately 2 kg
- Size : 257 (Width) x 371 (Depth) x 258 (Height) mm
- Material : Transparent panel : Antistatic plastic
Frame : Aluminum
- Applicable models (balances made by A&D)
 - BA-T/BA series (1 μ g model)
 - AD-4212A series
 - AD-4212B series
 - AD-4212C series
 - AD-4212D series

3. External dimensions

Caution

- Do not use organic solvents to clean the tabletop breeze break as that may damage the transparency and antistatic properties of the panels.
- Use a soft lint free cloth dampened with warm water and a mild detergent when cleaning the tabletop breeze break.