


HR-AZ/A・FX-CT シリーズ付属品変更のお知らせ

この度は、HR-AZ/A・FX-CTシリーズをお買い求め頂きありがとうございます。

HR-AZ/A・FX-CTシリーズの付属品及び風防内の寸法に関して、変更がありましたのでお知らせいたします。

1. 計量皿・皿受けが一個の計量皿ユニットになります。



計量皿・皿受け(別部品)


一個の計量皿ユニット(一部品)

変更前

変更後

取扱説明書 「製品構成(各部の名称)、設置と注意」 参照

2. 計量皿・皿受けが一個の計量皿ユニットになることによって、皿上高さと皿上内寸が変わります。




取扱説明書 「外形寸法図」 参照

Update to accessory for HR-AZ / A and FX-CT series



Thank you for purchasing the HR-AZ/A or FX-CT series. This is to notify customers of a change to an accessory for the series and to the dimensions inside the breeze break.

1. Both a weighing pan and a pan support are now provided together as a single weighing pan unit.



Weighing pan and pan support
(separate components)

Previous design




A single weighing pan unit that is composed
of the weighing pan and the pan support
*Please do not disassemble the unit.

Weighing pan unit
(one component)



New design

See 'Unpacking the Balance' and 'Precautions' in the instruction manual.

2. Accordingly, the height and inside dimension above the pan have been changed as shown below.



Previous dimensions



New dimensions

See 'External Dimensions' in the instruction manual.