

シリーズ名	天びんはかり		推奨する分銅 の商品コード	必要な数	参考	
	型名	商品コード			推奨分銅の質量	推奨分銅の精度等級
BA	BA-125D	BA125D	AD1603-100E2	1個	100g	E2
	BA-125DTE	BA125DTE	AD1603-100E2	1個	100g	E2
	BA-225	BA225	AD1603-200E2	1個	200g	E2
	BA-225D	BA225D	AD1603-200E2	1個	200g	E2
	BA-225DTE	BA225DTE	AD1603-200E2	1個	200g	E2
	BA-225TE	BA225TE	AD1603-200E2	1個	200g	E2
	BA-6DE	BA6DE	AD1603-5E2	1個	5g	E2
	BA-6DTE	BA6DTE	AD1603-5E2	1個	5g	E2
	BA-6E	BA6E	AD1603-5E2	1個	5g	E2
	BA-6TE	BA6TE	AD1603-5E2	1個	5g	E2
BM	BM-20	BM20	AD1603-20E2	1個	20g	E2
	BM-200	BM200	AD1603-200E2	1個	200g	E2
	BM-22	BM22	AD1603-20E2	1個	20g	E2
	BM-252	BM252	AD1603-200E2	1個	200g	E2
	BM-300	BM300	AD1603-200E2	1個	200g	E2
	BM-5	BM5	AD1603-5E2	1個	5g	E2
	BM-500	BM500	AD1603-200E2	1個	200g	E2
	BM-5D	BM5D	AD1603-5E2	1個	5g	E2
EJ-B(EJ)	EJ-1202B	EJ1202B-JA	AD1603-1KF1	1個	1000g	F1
	EJ-120B	EJ120B-JA	AD1603-100F2	1個	100g	F2
	EJ-123B	EJ123B-JA	AD1603-100F1	1個	100g	F1
	EJ-1500B	EJ1500B-JA	AD1603-1KF2	1個	1000g	F2
	EJ-2000B	EJ2000B-JA	AD1603-2KF2	1個	2000g	F2
	EJ-200B	EJ200B-JA	AD1603-200F2	1個	200g	F2
	EJ-3000B	EJ3000B-JA	AD1603-2KF2	1個	2000g	F2
	EJ-3002B	EJ3002B-JA	AD1603-2KF1	1個	2000g	F1
	EJ-300B	EJ300B-JA	AD1603-200F2	1個	200g	F2
	EJ-303B	EJ303B-JA	AD1603-200F1	1個	200g	F1
	EJ-4100B	EJ4100B-JA	AD1603-2KF1	2個	2000g	F1
	EJ-410B	EJ410B-JA	AD1603-200F2	2個	200g	F1
	EJ-54D2B	EJ54D2B-JA	AD1603-50F1	1個	50g	F1
	EJ-6100B	EJ6100B-JA	AD1603-5KF1	1個	5000g	F1
EJ-610B	EJ610B-JA	AD1603-500F1	1個	500g	F1	
EK-i	EK-1200i	EK1200i-JA	AD1603-1KF2	1個	1000g	F2
	EK-120i	EK120i-JA	AD1603-100F2	1個	100g	F2
	EK-12Ki	EK12Ki-JA	AD1603-10KF2	1個	10kg	F2
	EK-2000i	EK2000i-JA	AD1603-2KF2	1個	2000g	F2
	EK-200i	EK200i-JA	AD1603-200F2	1個	200g	F2
	EK-3000i	EK3000i-JA	AD1603-2KF2	1個	2000g	F2
	EK-300i	EK300i-JA	AD1603-200F2	1個	200g	F2
	EK-4100i	EK4100i-JA	AD1603-2KF1	2個	2000g	F1
	EK-410i	EK410i-JA	AD1603-200F1	2個	200g	F1
	EK-6000i	EK6000i-JA	AD1603-5KM1	1個	5000g	M1
	EK-600i	EK600i-JA	AD1603-500M1	1個	500g	M1
	EK-6100i	EK6100i-JA	AD1603-5KF1	1個	5000g	F1
	EK-610i	EK610i-JA	AD1603-500F1	1個	500g	F1
EK-L	EK-15KL	EK15KL-JA	AD1603-10KF1	1個	10kg	F1
	EK-30KL	EK30KL-JA	AD1603-20KF1	1個	20kg	F1
EW-i	EW-12Ki	EW12Ki-JA	AD1603-10KF2	1個	10kg	F2
	EW-1500i	EW1500i-JA	AD1603-1KF2	1個	1000g	F2
	EW-150i	EW150i-JA	AD1603-100F2	1個	100g	F2
FX-i	FX-1200i	FX1200i-JA	AD1603-1KF1	1個	1000g	F1
	FX-120i	FX120i-JA	AD1603-100F1	1個	100g	F1
	FX-2000i	FX2000i-JA	AD1603-2KF1	1個	2000g	F1
	FX-200i	FX200i-JA	AD1603-200F1	1個	200g	F1
	FX-3000i	FX3000i-JA	AD1603-2KF1	1個	2000g	F1
	FX-300i	FX300i-JA	AD1603-200F1	1個	200g	F1
	FX-5000i	FX5000i-JA	AD1603-5KF1	1個	5000g	F1
FX-500i	FX500i-JA	AD1603-500F1	1個	500g	F1	
FZ-i	FZ-1200i	FZ1200i-JA	AD1603-1KF1	1個	1000g	F1
	FZ-120i	FZ120i-JA	AD1603-100F1	1個	100g	F1
	FZ-2000i	FZ2000i-JA	AD1603-2KF1	1個	2000g	F1
	FZ-200i	FZ200i-JA	AD1603-200F1	1個	200g	F1
	FZ-3000i	FZ3000i-JA	AD1603-2KF1	1個	2000g	F1
	FZ-300i	FZ300i-JA	AD1603-200F1	1個	200g	F1
	FZ-5000i	FZ5000i-JA	AD1603-5KF1	1個	5000g	F1
FZ-500i	FZ500i-JA	AD1603-500F1	1個	500g	F1	

シリーズ名	天びんばかり		推奨する分銅 の商品コード	必要な数	参考	
	型名	商品コード			推奨分銅の質量	推奨分銅の精度等級
GF-A	GF-10001A	GF10001A	AD1603-10KF1	1個	10kg	F1
	GF-10002A	GF10002A	AD1603-10KF1	1個	10kg	F1
	GF-1003A	GF1003A	AD1603-1KF1	1個	1000g	F1
	GF-124A	GF124A	AD1603-100E2	1個	100g	E2
	GF-1603A	GF1603A	AD1603-1KF1	1個	1000g	F1
	GF-2002A	GF2002A	AD1603-2KF1	1個	2000g	F1
	GF-203A	GF203A	AD1603-200F1	1個	200g	F1
	GF-224A	GF224A	AD1603-200E2	1個	200g	E2
	GF-3002A	GF3002A	AD1603-2KF1	1個	2000g	F1
	GF-303A	GF303A	AD1603-200F1	1個	200g	F1
	GF-324A	GF324A	AD1603-200E2	1個	200g	E2
	GF-4002A	GF4002A	AD1603-2KF1	2個	2000g	F1
	GF-403A	GF403A	AD1603-200F1	2個	200g	F1
	GF-6001A	GF6001A	AD1603-5KF1	1個	5000g	F1
	GF-6002A	GF6002A	AD1603-5KF1	1個	5000g	F1
GF-603A	GF603A	AD1603-500F1	1個	500g	F1	
GH	GH-120	GH120	AD1603-100E2	1個	100g	E2
	GH-200	GH200	AD1603-200E2	1個	200g	E2
	GH-202	GH202	AD1603-200E2	1個	200g	E2
	GH-252	GH252	AD1603-200E2	1個	200g	E2
	GH-300	GH300	AD1603-200E2	1個	200g	E2
GP	GP-100K	GP100K	AD1603-20KF1	5個	20kg	F1
	GP-102K	GP102K	AD1603-20KF1	5個	20kg	F1
	GP-12K	GP12K	AD1603-10KF1	1個	10kg	F1
	GP-20K	GP20K	AD1603-20KF1	1個	20kg	F1
	GP-30K	GP30K	AD1603-20KF1	1個	20kg	F1
	GP-32K	GP32K	AD1603-20KF1	1個	20kg	F1
	GP-40K	GP40K	AD1603-20KF1	2個	20kg	F1
	GP-60K	GP60K	AD1603-20KF1	3個	20kg	F1
GP-61K	GP61K	AD1603-20KF1	3個	20kg	F1	
GR	GR-120	GR120	AD1603-100E2	1個	100g	E2
	GR-200	GR200	AD1603-200E2	1個	200g	E2
	GR-202	GR202	AD1603-200E2	1個	200g	E2
	GR-300	GR300	AD1603-200E2	1個	200g	E2
	GR-60	GR60	AD1603-50E2	1個	50g	E2
GX	GX-1000	GX1000	AD1603-1KF1	1個	1000g	F1
	GX-200	GX200	AD1603-200F1	1個	200g	F1
	GX-2000	GX2000	AD1603-2KF1	1個	2000g	F1
	GX-300	GX300	AD1603-200F1	1個	200g	F1
	GX-3000	GX3000	AD1603-2KF1	1個	2000g	F1
	GX-400	GX400	AD1603-200F1	2個	200g	F1
	GX-4000	GX4000	AD1603-2KF1	2個	2000g	F1
	GX-600	GX600	AD1603-500F1	1個	500g	F1
	GX-6000	GX6000	AD1603-5KF1	1個	5000g	F1
	GX-6100	GX6100	AD1603-5KF1	1個	5000g	F1
	GX-800	GX800	AD1603-500F1	1個	500g	F1
GX-8000	GX8000	AD1603-5KF1	1個	5000g	F1	
GX-AE/A	GX-10001A	GX10001A	AD1603-10KF1	1個	10kg	F1
	GX-10002A	GX10002A	AD1603-10KF1	1個	10kg	F1
	GX-1003A	GX1003A	AD1603-1KF1	1個	1000g	F1
	GX-124A	GX124A	AD1603-100E2	1個	100g	E2
	GX-1603A	GX1603A	AD1603-1KF1	1個	1000g	F1
	GX-2002A	GX2002A	AD1603-2KF1	1個	2000g	F1
	GX-203A	GX203A	AD1603-200F1	1個	200g	F1
	GX-224A	GX224A	AD1603-200E2	1個	200g	E2
	GX-3002A	GX3002A	AD1603-2KF1	1個	2000g	F1
	GX-303A	GX303A	AD1603-200F1	1個	200g	F1
	GX-324A	GX324A	AD1603-200E2	1個	200g	E2
	GX-4002A	GX4002A	AD1603-2KF1	2個	2000g	F1
	GX-403A	GX403A	AD1603-200F1	2個	200g	F1
	GX-6001A	GX6001A	AD1603-5KF1	1個	5000g	F1
	GX-6002A	GX6002A	AD1603-5KF1	1個	5000g	F1
	GX-603A	GX603A	AD1603-500F1	1個	500g	F1
	GX-124AE	GX124AE	AD1603-100E2	1個	100g	E2
	GX-224AE	GX224AE	AD1603-200E2	1個	200g	E2
GX-324AE	GX324AE	AD1603-200E2	1個	200g	E2	

シリーズ名	天びんばかり		推奨する分銅 の商品コード	必要な数	参考	
	型名	商品コード			推奨分銅の質量	推奨分銅の精度等級
GX-K	GX-10K	GX10K	AD1603-10KF1	1個	10kg	F1
	GX-12K	GX12K	AD1603-10KF1	1個	10kg	F1
	GX-20K	GX20K	AD1603-20KF1	1個	20kg	F1
	GX-30K	GX30K	AD1603-20KF1	1個	20kg	F1
	GX-32K	GX32K	AD1603-20KF1	1個	20kg	F1
	GX-8K	GX8K	AD1603-5KF1	1個	5000g	F1
	GX-8K2	GX8K2	AD1603-5KF1	1個	5000g	F1
GF-M	GF-10202M	GF10202M	AD1603-10KF1	1個	10kg	F1
	GF-12001M	GF12001M	AD1603-10KF1	1個	10kg	F1
	GF-22001M	GF22001M	AD1603-20KF1	1個	20kg	F1
	GF-32001M	GF32001M	AD1603-20KF1	1個	20kg	F1
	GF-32001MD	GF32001MD	AD1603-20KF1	1個	20kg	F1
	GF-8202M	GF8202M	AD1603-5KF1	1個	5000g	F1
	GF-8202MD	GF8202MD	AD1603-5KF1	1個	5000g	F1
GX-M	GX-10202M	GX10202M	AD1603-10KF1	1個	10kg	F1
	GX-12001M	GX12001M	AD1603-10KF1	1個	10kg	F1
	GX-22001M	GX22001M	AD1603-20KF1	1個	20kg	F1
	GX-32001M	GX32001M	AD1603-20KF1	1個	20kg	F1
	GX-32001MD	GX32001MD	AD1603-20KF1	1個	20kg	F1
	GX-8202M	GX8202M	AD1603-5KF1	1個	5000g	F1
	GX-8202MD	GX8202MD	AD1603-5KF1	1個	5000g	F1
HR/HR-A	HR-100A(AZ)	HR100AZ-JA	AD1603-100E2	1個	100g	E2
	HR-120	HR120	AD1603-100E2	1個	100g	E2
	HR-150A(AZ)	HR150AZ-JA	AD1603-100E2	1個	100g	E2
	HR-200	HR200	AD1603-200E2	1個	200g	E2
	HR-202i	HR202i	AD1603-200E2	1個	200g	E2
	HR-250A(AZ)	HR250AZ-JA	AD1603-200E2	1個	200g	E2
	HR-251A(AZ)	HR251AZ-JA	AD1603-200E2	1個	200g	F1
	HR-300i	HR300i	AD1603-200E2	1個	200g	E2
HR-60	HR60	AD1603-50E2	1個	50g	E2	
MC	MC-1000	MC1000	AD1603-1KE2	1個	1000g	E2
	MC-100KS	MC100KS	AD1603-20KF1	5個	20kg	F1
	MC-10K	MC10K	AD1603-10KF1	1個	10kg	F1
	MC-30K	MC30K	AD1603-20KF1	1個	20kg	F1
	MC-6100	MC6100	AD1603-5KF1	1個	5000g	F1